

Resuscitation sheet

in keychain




Well protected in an emergency

The first aider's self-protection has the highest priority in an emergency. With our HERZmed resuscitation sheet, you can practice life-saving first aid measures safely and without endangering yourself in an emergency. Stowed in a robust nylon protective bag with keyring, the sheet can be conveniently and discreetly attached to a keychain or backpack and is always at hand. The protective cloth is fixed over the patient's face during ventilation measures in the context of cardiopulmonary resuscitation and thus forms a protective barrier between the patient and the rescuer, which reduces the contact inhibition threshold. The valve incorporated into the sheet minimizes the risk of germ transmission.

An instruction manual printed on the sheet guides the user through the individual steps of the ventilation measures in an easy-to-understand manner.

Details



Valve



Keyring



Printed instruction



Highlights

- Safe protection against infection transmissions
- Reduction of inhibition threshold through physical contact
- Always ready to hand on keyring, backpack, bag & Co.
- Handy and small
- With easy to understand, printed operating instruction
- Can be customized with company logo or HERZmed logo
- Bag available in eight different, exciting colors

Latex-free

The ventilation sling is latex-free. You can find the exact materials in the technical specifications.

Suitable for advertising

With this smart advertising medium you are sure to be successful and do good on top of it. Already from a quantity of 500 pieces.

Velcro fastener

The Velcro ensures secure closure and quick opening in case of emergency.

Instruction included

The instructions are printed directly on the ventilation sheet to make the application as efficient as possible.

Robust

The protective bag is made of sturdy nylon material and protects the respirator sheet from external damage.

Always present

Thanks to the handy design and key ring, the ventilation sheet is always with you and ready for use.



In Detail

The keychain offers the possibility in an emergency during mouth-to-mouth resuscitation of the self-protection of the first secure helpers. Optimally packed in a nylon protective bag, this little lifesaver can be taken anywhere. An integrated, flat valve offers the first aider safe hygiene protection. Easy-to-understand instructions for cardiopulmonary resuscitation with illustrations can be found on the front of the ventilation foil.

Technical specifications

Size:

Schutztasche: 50 x 50 x 10 mm (L x B x H)

Tuch, offen: 280 x 200 (L x B)

Tuch, geschlossen: 50 x 50 x 10 mm (L x B x H)

Weight:

Ca. 10 g

Coloring:

See below

Packaging:

Individually packed in protective bag

Material:

Beatmungstuch: Polyvinylchlorid (PVC)

Ventil: Polyvinylchlorid (PVC)

Filtermedium: Hydrophobe Membran

Sterility:

No

Order information

Customized logo

Art.-No.	Description	Image
HZ-6004	Keychain green, customized logo	
HZ-6005	Keychain light pink, customized logo	
HZ-6006	Keychain yellow, customized logo	
HZ-6007	Keychain red, customized logo	
HZ-6008	Keychain white, customized logo	
HZ-6009	Keychain purple, customized logo	
HZ-6010	Keychain black, customized logo	
HZ-6011	Keychain blue, customized logo	
HZ-6023	Keychain orange, customized logo	

Minimum order quantity 1 PU = 500 pieces:

Customized printed nylon bag of the respiratory wipes in the keychain possible from a PU of 500 pieces.

HERZmed logo

Art.-No.	Description	Image
HZ-6015	Keychain blue, HERZmed logo	
HZ-6016	Keychain red, HERZmed logo	
HZ-6017	Keychain black, HERZmed logo	

No minimum order required.

Order and advice by:



HERZmed GmbH

Herkulan-Schwaiger-Gasse 16
82487 Oberammergau

Tel.: +49 (0) 911 / 14 88 959 - 0
Fax: +49 (0) 911 / 14 88 959 - 9

Mail: herzmed@herzmed.de
Web: www.herzmed.de

The products shown are only an excerpt of our product range. Exact details and product information are available on request or at www.herzmed.de. All prices marked with an asterisk are exclusive of VAT. All brand names used remain the property of the respective brand owner. We assume no liability for errors, possible misprints and color deviations. Our general terms and conditions apply, which you can also view at www.herzmed.de.

