

# Rescue sheet

*ONE-WAY BLACK*



## Carries more than responsibility

The HERZmed EINWEG rescue sling is suitable for transporting people from rough terrain or difficult to access areas. It can be used, for example, by rescue services, the fire department or mountain rescue services.

Load straps, which are sewn crosswise to the sling, provide high stability. The HERZmed EINWEG rescue sling is used for transporting patients lying down.

It enables patients to be transported safely and comfortably as if wrapped in a carry bag. This type of sheet is particularly useful in emergency situations as it is hygienic and easy to handle. It can be easily disposed of after use.

The HERZmed EINWEG rescue sheet is made of high-strength polyester fabric coated on both sides with PVC. The load straps are made of woven polyester covered with PE plastic tubing.

## Details



**Additional protection with headboard**



**and footmuff**



**double PVC-coated polyester tarpaulin material**

## Highlights

- With leg bag and head protection
- Water-resistant and hard-wearing PVC coated polyester tarpaulin material
- 10/14 handles
- Size of the lying surface: L 200 cm x 80 cm
- Size of the lying surface: XL 200 cm x 98 cm
- Load capacity L 200 kg or XL 250 kg
- Rescue sheet complies with DIN EN 1865

## Lightweight and compact

Space-saving and lightweight.

## Headboard

Integrated headboard for protection.

## High load capacity

Load capacity of up to  
L = 200 kg  
XL = 250 kg

## Size L 10 / XL14 Straps

The straps are made of woven polyester material covered with PE plastic tubing.

## Large lying surface

Version L = 200 x 80 cm  
Version XL = 200 x 98 cm.

## Durable

Made from hard-wearing coated on both sides with PVC coated high-strength polyester fabric.

Impermeable to liquids and 100% X-ray transparent.

## Leg bag

The leg bag prevents the patient from slipping out.

## Versatile in use

Load straps, which are sewn crosswise to the sling, give the HERZmed rescue sling high stability. The HERZmed EINWEG rescue sling is suitable for transporting people from rough terrain or areas that are difficult to access. It can be used, for example, by rescue services, the fire department or mountain rescue services. Dispose of in normal household waste as residual waste.

### Available in:

#### Size L



#### Size XL



## Technical specifications

### Size:

Version L = 200 x 80 cm  
Version XL = 200 x 98 cm

### Load capacity:

Version L = 200 kg  
Version XL = 250 kg

**Colour:** BLACK

### Packaging:

**Individually packed in transparent foil**

**Sterile:** No

**Storage:** Store in a dry place

**Disposal:** In normal household waste as residual waste

### Materials:

Lying surface:  
Made of hard-wearing high-strength polyester fabric coated on both sides with PVC. Impermeable to liquids and 100% X-ray transparent.

Carrying straps made of durable woven polyester covered with PE plastic tubing.

## Order information:

Art.-No.	Description	Image
190458	Rescue sheets ONE-WAY BLACK L	
190459	Rescue sheets ONE-WAY BLACK XL	

## Order and advice:

**Vertrieb/Sales Office:**  
HERZmed GmbH  
Rottenbacher Straße 9  
82497 Unterammergau

**Tel.:** +49 (0) 8822 824 0002  
**Fax:** +49 (0) 8822 824 0006

**Verwaltung/Management:**  
HERZmed GmbH  
Herkulan-Schwaiger-Gasse 16  
82487 Oberammergau

**Mail:** herzmed@herzmed.de  
**Web:** www.herzmed.com

The products shown are only an excerpt of our product range. Exact details and product information are available on request or at [www.herzmed.de](http://www.herzmed.de). All prices marked with an asterisk are exclusive of VAT. All brand names used remain the property of the respective brand owner. We assume no liability for errors, possible misprints and color deviations. Our general terms and conditions apply, which you can also view at [www.herzmed.de](http://www.herzmed.de).

